

EMS hotplate programming instructions

To create temperature setpoints

Push P. Key for about 5 secs. (Green display will show decimal point)

Push P Key again. (Green display shows 0)

Push up key until green display reads 381

Push P (you are now into the program menu)

Cycle through program menu using the up key

NSP (number of set points) Change to required no.

Push P now set the temperatures required

Then set temperatures of each set point

All you need to do now is go to PAN menu and change USrb to CHSP

(Change set point using the U Key)

By holding down the up key after you have changed all your settings. It will send you back to the run conditions.

To adjust setpoint temperature

Press P

Use up/down arrows to adjust to new temperature

Press P again

Plate will now ramp to new temperature\